




NXPE2 mouse mAb

Catalog No	YP-mAb-11677
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NXPE2 FAM55B
Immunogen	Synthesized peptide derived from human NXPE2 AA range: 8-58
Specificity	This antibody detects endogenous levels of NXPE2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM55B;FAM55B Ens;NXPE2;NXPE2 FAM55B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., similarity:Belongs to the FAM55 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NXPE2 mouse mAb