




CK048 mouse mAb

Catalog No	YP-mAb-11690
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C11orf48
Immunogen	Synthesized peptide derived from human CK048. AA range 37-67
Specificity	This antibody detects endogenous levels of CK048 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C11orf48;CK048
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in bladder cancer tissues (at protein level).
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CK048 mouse mAb