




SBNO1 mouse mAb

Catalog No	YP-mAb-11701
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SBNO1 MOP3
Immunogen	Synthesized peptide derived from human SBNO1 AA range: 383-433
Specificity	This antibody detects endogenous levels of SBNO1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Monocyte protein 3;MOP 3;MOP3;SBNO1;Sno;SBNO1 MOP3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the SBNO family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SBNO1 mouse mAb