




FA63B mouse mAb

Catalog No	YP-mAb-11713
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM63B KIAA1164
Immunogen	Synthesized peptide derived from human FA63B AA range: 302-352
Specificity	This antibody detects endogenous levels of FA63B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA63B;FAM63B KIAA1164
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleoplasm,
Tissue Specificity	
Function	similarity:Belongs to the FAM63 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA63B mouse mAb