




S35F3 mouse mAb

Catalog No	YP-mAb-11734
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC35F3
Immunogen	Synthesized peptide derived from human S35F3 AA range: 305-355
Specificity	This antibody detects endogenous levels of S35F3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S35F3;SLC35F3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	function:Putative solute transporter .,similarity:Belongs to the SLC35F solute transporter family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using S35F3 mouse mAb