




FA45A mouse mAb

Catalog No	YP-mAb-11736
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM45A
Immunogen	Synthesized peptide derived from human FA45A AA range: 126-176
Specificity	This antibody detects endogenous levels of FA45A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA45A;FAM45A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Late endosome .
Tissue Specificity	
Function	similarity:Belongs to the FAM45 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA45A mouse mAb