




LACB2 mouse mAb

Catalog No	YP-mAb-11750
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LACTB2 CGI-83
Immunogen	Synthesized peptide derived from human LACB2 AA range: 19-69
Specificity	This antibody detects endogenous levels of LACB2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.1.27.-;CGI-83;CGI83;CGI 83;Beta-lactamase-like protein 2;LACB2;LACTB2 CGI-83
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	similarity:Belongs to the metallo-beta-lactamase superfamily. Glyoxalase II family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LACB2 mouse mAb