




TM86B mouse mAb

Catalog No	YP-mAb-11793
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM86B
Immunogen	Synthesized peptide derived from human TM86B AA range: 46-96
Specificity	This antibody detects endogenous levels of TM86B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM86B;TMEM86B
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Membrane ; Multi-pass membrane protein . Cytoplasm .
Tissue Specificity	
Function	similarity:Belongs to the TMEM86 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TM86B mouse mAb