




TTC12 mouse mAb

Catalog No	YP-mAb-11801
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC12
Immunogen	Synthesized peptide derived from human TTC12 AA range: 579-629
Specificity	This antibody detects endogenous levels of TTC12 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ13859;FLJ20535;TPARM;TPR repeat protein 12;TTC12
Observed Band	78kD
Calculated Molecular Weight	78kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in testis and in epithelial cells of trachea and bronchial tube.
Function	similarity:Contains 3 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4°over night