




GATL2 mouse mAb

Catalog No	YP-mAb-11803
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GATSL2
Immunogen	Synthesized peptide derived from human GATL2 AA range: 79-129
Specificity	This antibody detects endogenous levels of GATL2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GATS like protein 2;GATSL2;GATL2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GATL2 mouse mAb