




RCOR3 mouse mAb

Catalog No	YP-mAb-11810
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RCOR3 KIAA1343
Immunogen	Synthesized peptide derived from human RCOR3 AA range: 318-368
Specificity	This antibody detects endogenous levels of RCOR3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1343;RCOR3;REST corepressor 3;RP11 318L16.1;RCOR3 KIAA1343
Observed Band	
Calculated Molecular Weight	54kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May act as a component of a corepressor complex that represses transcription .,similarity:Belongs to the CoREST family.,similarity:Contains 1 ELM2 domain.,similarity:Contains 2 SANT domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RCOR3 mouse mAb