




GG6L6 mouse mAb

Catalog No	YP-mAb-11816
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GOLGA6L6
Immunogen	Synthesized peptide derived from human GG6L6 AA range: 26-76
Specificity	This antibody detects endogenous levels of GG6L6 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GG6L6;GOLGA6L6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,similarity:Belongs to the GOLGA6 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GG6L6 mouse mAb