




PAQRA mouse mAb

Catalog No	YP-mAb-11823
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MMD2 PAQR10
Immunogen	Synthesized peptide derived from human PAQRA AA range: 30-80
Specificity	This antibody detects endogenous levels of PAQRA at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAQRA;MMD2 PAQR10
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Shows restricted expression with highest levels in brain and testis.
Function	similarity:Belongs to the ADIPOR family.,
Background	This gene encodes a member of the PAQR (progesterone and adipoQ receptor) family. Members of this family are evolutionarily conserved with significant sequence identity to bacterial hemolysin-like proteins and are defined by a set of seven transmembrane domains. The protein encoded by this gene localizes to the Golgi apparatus to modulate Ras signaling. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jun 2012],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using PAQRA mouse mAb