




SPI2B mouse mAb

Catalog No	YP-mAb-11846
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPIN2B SPIN2
Immunogen	Synthesized peptide derived from human SPI2B AA range: 73-123
Specificity	This antibody detects endogenous levels of SPI2B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPIN-2B;SPIN2_duplicate;SPIN-2;SPIN2;SPIN 2B;SPI2B;SPIN2B SPIN2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Detected in all the examined tissues with highest expression in liver, followed by heart, stomach, kidney, skeletal muscle, placenta, and pancreas.
Function	function:Involved in the regulation of cell cycle progression, this activity is related to the inhibition of apoptosis following the removal of essential growth factors., miscellaneous:Overexpression in murine myeloid cell line 32Dcl3 causes G2/M arrest.,similarity:Belongs to the SPIN/STSY family.,tissue specificity:Detected in all the examined tissues with highest expression in liver, followed by heart, stomach, kidney, skeletal muscle, placenta, and pancreas.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SPI2B mouse mAb