




AN13D mouse mAb

Catalog No	YP-mAb-11868
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANKRD13D
Immunogen	Synthesized peptide derived from human AN13D AA range: 382-432
Specificity	This antibody detects endogenous levels of AN13D at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKRD13D;AN13D
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane . Late endosome . Interaction with EGFR may enhance association with the cell membrane.
Tissue Specificity	
Function	similarity:Contains 3 UIM (ubiquitin-interacting motif) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using AN13D mouse mAb