




VGLL1 mouse mAb

Catalog No	YP-mAb-11877
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VGLL1 TDU
Immunogen	Synthesized peptide derived from human VGLL1 AA range: 84-134
Specificity	This antibody detects endogenous levels of VGLL1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein TONDU;TDU;vestigial like 1 (Drosophila);Vgl 1;VGL1;VGLL1;VGLL1 TDU
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	function:May act as a specific coactivator for the mammalian TEFs., similarity:Belongs to the vestigial family.,subunit:Interacts with TEFs.,
Background	The protein encoded by this gene binds proteins of the TEA domain family of transcription factors (TEFs) through the Vg (vestigial) homology region found in its N-terminus. It may thus function as a specific coactivator for the mammalian TEFs. [provided by RefSeq, Sep 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using VGLL1 mouse mAb