




SYPL2 mouse mAb

Catalog No	YP-mAb-11935
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SYPL2
Immunogen	Synthesized peptide derived from human SYPL2 AA range: 160-210
Specificity	This antibody detects endogenous levels of SYPL2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MG29;mitsugumin 29;synaptophysin like 2;Synaptophysin-like protein 2;SYPL2; SYPL2_HUMAN
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Triad junction, the junctional complex between the transverse tubule and the sarcoplasmic reticulum. .
Tissue Specificity	
Function	function:Involved in communication between the T-tubular and junctional sarcoplasmic reticulum (SR) membranes.,similarity:Belongs to the synaptophysin/synaptobrevin family.,similarity:Contains 1 MARVEL domain., subcellular location:Triad junction, the junctional complex between the transverse tubule and the sarcoplasmic reticulum.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SYPL2 mouse mAb