




RTKN2 mouse mAb

Catalog No	YP-mAb-11944
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RTKN2 PLEKHK1
Immunogen	Synthesized peptide derived from human RTKN2 AA range: 473-523
Specificity	This antibody detects endogenous levels of RTKN2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PH domain-containing family K member 1;Pleckstrin homology domain-containing family K;RTKN2;RTKN2 PLEKHK1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.
Function	function:May play an important role in lymphopoiesis.,induction:Up-regulated in cells resistant to 25-hydroxy cholesterol (25-OHC). Down-regulated in lymphocytes activated by treatment with phytohemagglutinin (PHA)., similarity:Contains 1 PH domain.,tissue specificity:Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RTKN2 mouse mAb