




KDEL2 mouse mAb

Catalog No	YP-mAb-11971
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KDELC2 UNQ1904/PRO4350
Immunogen	Synthesized peptide derived from human KDEL2 AA range: 414-464
Specificity	This antibody detects endogenous levels of KDEL2 at Human/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein O-xylosyltransferase POGLUT3;Protein O-glucosyltransferase 3; POGLUT3;KDEL motif-containing protein 2;EC:2.4.2.-;KDEL2;KDELC2 UNQ1904/PRO4350
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum lumen .
Tissue Specificity	
Function	similarity:Belongs to the KDELC family.,similarity:Contains 1 filamin repeat.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KDEL2 mouse mAb