




ANGL6 mouse mAb

Catalog No	YP-mAb-11981
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANGPTL6 AGF ARP5 UNQ152/PRO178
Immunogen	Synthesized peptide derived from human ANGL6 AA range: 216-266
Specificity	This antibody detects endogenous levels of ANGL6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANGL6;ANGPTL6 AGF ARP5 UNQ152/PRO178
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:May play a role in the wound healing process. May promote epidermal proliferation, remodeling and regeneration. May promote the chemotactic activity of endothelial cells and induce neovascularization. May counteract high-fat diet-induced obesity and related insulin resistance through increased energy expenditure.,similarity:Contains 1 fibrinogen C-terminal domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ANGL6 mouse mAb