



IAH1 mouse mAb

Catalog No	YP-mAb-11984
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IAH1
Immunogen	Synthesized peptide derived from human IAH1 AA range: 22-72
Specificity	This antibody detects endogenous levels of IAH1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IAH1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular exosome,
Tissue Specificity	
Function	function:Probable lipase.,similarity:Belongs to the 'GDSL' lipolytic enzyme family. IAH1 subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using IAH1 mouse mAb

