




SHF mouse mAb

Catalog No	YP-mAb-11986
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHF
Immunogen	Synthesized peptide derived from human SHF AA range: 163-213
Specificity	This antibody detects endogenous levels of SHF at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SHF
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in skeletal muscle, brain, liver, prostate, testis, ovary, small intestine and colon.
Function	function:Adapter protein which may play a role in the regulation of apoptosis in response to PDGF.,PTM:May become phosphorylated upon binding to PDGFRA.,similarity:Contains 1 SH2 domain.,subunit:Interacts with phosphorylated 'Tyr-720' of PDGFRA via its SH2 domain.,tissue specificity:Expressed in skeletal muscle, brain, liver, prostate, testis, ovary, small intestine and colon.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SHF mouse mAb