




2B1G mouse mAb

Catalog No	YP-mAb-11995
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HLA-DRB1
Immunogen	Synthesized peptide derived from human 2B1G AA range: 159-209
Specificity	This antibody detects endogenous levels of 2B1G at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Human leukocyte antigen DRB1;HLA DRB1;HLA DRB;HLA DR1B;HLA class II histocompatibility antigen;DRB1 beta chain;2B1G;HLA-DRB1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using 2B1G mouse mAb