




TMTC4 mouse mAb

Catalog No	YP-mAb-12000
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMTC4
Immunogen	Synthesized peptide derived from human TMTC4 AA range: 28-78
Specificity	This antibody detects endogenous levels of TMTC4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMTC4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Endoplasmic reticulum .
Tissue Specificity	
Function	similarity:Belongs to the TMTC family.,similarity:Contains 8 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TMTC4 mouse mAb