




## TVC mouse mAb

|                                    |                                                                                                                                                                                                           |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-12005                                                                                                                                                                                              |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                       |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                           |
| <b>Applications</b>                | WB                                                                                                                                                                                                        |
| <b>Gene Name</b>                   | TRGV3 TCRB                                                                                                                                                                                                |
| <b>Immunogen</b>                   | Synthesized peptide derived from human TVC AA range: 79-129                                                                                                                                               |
| <b>Specificity</b>                 | This antibody detects endogenous levels of TVC at Human                                                                                                                                                   |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                   |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                     |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                      |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                           |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                   |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                      |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                              |
| <b>Synonyms</b>                    | TVC;TRGV3 TCRB                                                                                                                                                                                            |
| <b>Observed Band</b>               |                                                                                                                                                                                                           |
| <b>Calculated Molecular Weight</b> |                                                                                                                                                                                                           |
| <b>Cell Pathway</b>                | Cell membrane .                                                                                                                                                                                           |
| <b>Tissue Specificity</b>          |                                                                                                                                                                                                           |
| <b>Function</b>                    | miscellaneous:This is an aberrantly rearranged nucleotide sequence. Two bases are deleted at nucleotide position 357 in this clone. However, the protein sequence can be deduced from the mRNA sequence., |
| <b>Background</b>                  |                                                                                                                                                                                                           |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                                                      |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                         |



## Products Images



Western Blot analysis of various cells using TVC mouse mAb