




BIN3 mouse mAb

Catalog No	YP-mAb-12012
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BIN3
Immunogen	Synthesized peptide derived from human BIN3 AA range: 184-234
Specificity	This antibody detects endogenous levels of BIN3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bridging integrator 3;BIN3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	Ubiquitously expressed except in brain.
Function	function:Involved in cytokinesis and septation where it has a role in the localization of F-actin.,similarity:Contains 1 BAR domain.,tissue specificity:Ubiquitously expressed except in brain.,
Background	The product of this gene is a member of the BAR domain protein family. The encoded protein is comprised solely of a BAR domain which is predicted to form coiled-coil structures and proposed to mediate dimerization, sense and induce membrane curvature, and bind small GTPases. BAR domain proteins have been implicated in endocytosis, intracellular transport, and a diverse set of other processes. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using BIN3 mouse mAb