




FKB15 mouse mAb

Catalog No	YP-mAb-12024
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FKBP15 KIAA0674
Immunogen	Synthesized peptide derived from human FKB15 AA range: 847-897
Specificity	This antibody detects endogenous levels of FKB15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	133 kDa FK506 binding protein;133 kDa FKBP;FK506 binding protein 15;FKBP 133;FKBP 15;FKB15;FKBP15 KIAA0674
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cell projection, axon . Early endosome . Present in axons and neuronal growth cones. .
Tissue Specificity	
Function	domain:The PPlase FKBP-type domain seems to be inactive both for FK506-binding and enzymatic activity.,function:May be involved in the cytoskeletal organization of neuronal growth cones. Seems to be inactive as a PPlase.,similarity:Belongs to the FKBP-type PPlase family.,similarity:Contains 1 PPlase FKBP-type domain.,subcellular location:Present in axons and neuronal growth cones.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FKB15 mouse mAb