




KRBA1 mouse mAb

Catalog No	YP-mAb-12027
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KRBA1 KIAA1862
Immunogen	Synthesized peptide derived from human KRBA1 AA range: 738-788
Specificity	This antibody detects endogenous levels of KRBA1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KRBA1;KRBA1 KIAA1862
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,
Tissue Specificity	Expressed in brain (cerebellum).
Function	tissue specificity:Expressed in brain (cerebellum).,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KRBA1 mouse mAb