




PCGF5 mouse mAb

Catalog No	YP-mAb-12033
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PCGF5 RNF159
Immunogen	Synthesized peptide derived from human PCGF5 AA range: 65-115
Specificity	This antibody detects endogenous levels of PCGF5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCGF5;PCGF5 RNF159
Observed Band	
Calculated Molecular Weight	28kD
Cell Pathway	Nucleus . Nucleus, nucleoplasm . Recruited by the non-coding RNA Xist to specific nuclear foci that probably correspond to the inactivated X chromosome. .
Tissue Specificity	
Function	function:Probable component of some Polycomb group (PcG) multiprotein complex, a complex required to maintain the transcriptionally repressive state of some genes.,similarity:Contains 1 RING-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PCGF5 mouse mAb