




DJC21 mouse mAb

Catalog No	YP-mAb-12045
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC21 DNAJA5
Immunogen	Synthesized peptide derived from human DJC21 AA range: 255-305
Specificity	This antibody detects endogenous levels of DJC21 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein GS3;GS3;DNAJA5;DJC21;DNAJC21 DNAJA5
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Cytoplasm . Nucleus . Nucleus, nucleolus . Within the nucleus, localizes primarily to the nucleolus. .
Tissue Specificity	Expressed in brain, placenta, kidney and pancreas.
Function	function:May act as a co-chaperone for HSP70.,similarity:Contains 1 J domain., similarity:Contains 2 C2H2-type zinc fingers.,tissue specificity:Brain, placenta, kidney and pancreas.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night