




POTEC mouse mAb

Catalog No	YP-mAb-12051
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POTEC A26B2 POTE18
Immunogen	Synthesized peptide derived from human POTEC AA range: 265-315
Specificity	This antibody detects endogenous levels of POTEC at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POTEC;POTEC A26B2 POTE18
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in prostate and testis.
Function	similarity:Belongs to the POTE family.,similarity:Contains 6 ANK repeats.,similarity:Contains 7 ANK repeats.,tissue specificity:Expressed in prostate and testis.,tissue specificity:Expressed in prostate, ovary, testis, placenta and prostate cancer cell lines. Localizes to basal and terminal prostate epithelial cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using POTEC mouse mAb