




ZIM3 mouse mAb

Catalog No	YP-mAb-12095
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZIM3 ZNF657
Immunogen	Synthesized peptide derived from human ZIM3 AA range: 98-148
Specificity	This antibody detects endogenous levels of ZIM3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZIM 3;Zinc finger protein 657;ZNF657;ZIM3;ZIM3 ZNF657
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 11 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZIM3 mouse mAb