




PALMD mouse mAb

Catalog No	YP-mAb-12118
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PALMD C1orf11 PALML
Immunogen	Synthesized peptide derived from human PALMD AA range: 4-54
Specificity	This antibody detects endogenous levels of PALMD at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf11;Palmd;PALMD_HUMAN;Palmdelphin;PALML;Paralemmin-Like Protein; Paralemnin-like;Paralemnin-like protein;PALMD C1orf11 PALML
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cell projection, dendrite . Cell projection, dendritic spine .
Tissue Specificity	Ubiquitous. Most abundant in cardiac and skeletal muscle.
Function	sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the paralemmin family.,tissue specificity:Ubiquitous. Most abundant in cardiac and skeletal muscle.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PALMD mouse mAb