



TIGD7 mouse mAb


Catalog No	YP-mAb-12129
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TIGD7
Immunogen	Synthesized peptide derived from human TIGD7 AA range: 308-358
Specificity	This antibody detects endogenous levels of TIGD7 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIGD 7;Sancho;TIGD7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in all tissues tested. Higher expression in testis and ovary.
Function	similarity:Belongs to the tigger transposable element derived protein family., similarity:Contains 1 DDE domain.,similarity:Contains 1 HTH CENPB-type DNA-binding domain.,similarity:Contains 1 HTH psq-type DNA-binding domain.,tissue specificity:Expressed in all tissues tested. Higher expression in testis and ovary.,
Background	The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TIGD7 mouse mAb