




GGNB2 mouse mAb

Catalog No	YP-mAb-12151
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GGNBP2 LCRG1 LZK1 ZNF403
Immunogen	Synthesized peptide derived from human GGNB2 AA range: 225-275
Specificity	This antibody detects endogenous levels of GGNB2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GGNB2;GGNBP2 LCRG1 LZK1 ZNF403
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle . Associated with vesicular structures. .
Tissue Specificity	Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Expressed more abundantly in heart, pancreas and skeletal muscle.
Function	function:May be involved in spermatogenesis.,miscellaneous:Strongly down-regulated in 40% of primary laryngeal carcinoma and in 6 of 10 various cancer cell lines.,subcellular location:Associated with vesicular structures., subunit:Interacts with GGN.,tissue specificity:Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Expressed more abundantly in heart, pancreas and skeletal muscle.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GGNB2 mouse mAb