




SOLH1 mouse mAb

Catalog No	YP-mAb-12178
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SOHLH1 C9orf157 NOHLH TEB2
Immunogen	Synthesized peptide derived from human SOLH1 AA range: 94-144
Specificity	This antibody detects endogenous levels of SOLH1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf157;NOHLH;SOHLH1;TEB2;SOLH1;SOHLH1 C9orf157 NOHLH TEB2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	function:Probable transcription factor required during spermatogenesis and oogenesis.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
Background	This gene encodes one of testis-specific transcription factors which are essential for spermatogenesis, oogenesis and folliculogenesis. This gene is located on chromosome 9. Mutations in this gene are associated with nonobstructive azoospermia. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SOLH1 mouse mAb