




REXO4 mouse mAb

Catalog No	YP-mAb-12194
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	REXO4 PMC2 XPMC2H
Immunogen	Synthesized peptide derived from human REXO4 AA range: 312-362
Specificity	This antibody detects endogenous levels of REXO4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	XPMC2;REX4;REX 4;PMC2;hPMC2;REXO4;REXO4 PMC2 XPMC2H
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	similarity:Belongs to the REXO4 family.,similarity:Contains 1 exonuclease domain.,subunit:Can bind ESR1 and ESR2. This interaction is abrogated by estrogen and augmented by tamoxifen treatment.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using REXO4 mouse mAb