




AIM1L mouse mAb

Catalog No	YP-mAb-12200
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AIM1L CRYBG2
Immunogen	Synthesized peptide derived from human AIM1L AA range: 212-262
Specificity	This antibody detects endogenous levels of AIM1L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AIM1L;AIM1L CRYBG2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 ricin B-type lectin domain.,similarity:Contains 10 beta/gamma crystallin 'Greek key' domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using AIM1L mouse mAb