




## ZN547 mouse mAb

|                             |                                                                                                                                                                                                                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-mAb-12212                                                                                                                                                                                                                                                                           |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                    |
| Reactivity                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                        |
| Applications                | WB                                                                                                                                                                                                                                                                                     |
| Gene Name                   | ZNF547                                                                                                                                                                                                                                                                                 |
| Immunogen                   | Synthesized peptide derived from human ZN547 AA range: 312-362                                                                                                                                                                                                                         |
| Specificity                 | This antibody detects endogenous levels of ZN547 at Human                                                                                                                                                                                                                              |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                |
| Source                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                  |
| Purification                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                   |
| Dilution                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                        |
| Concentration               | 1 mg/ml                                                                                                                                                                                                                                                                                |
| Purity                      | ≥90%                                                                                                                                                                                                                                                                                   |
| Storage Stability           | -20°C/1 year                                                                                                                                                                                                                                                                           |
| Synonyms                    | ZN547;ZNF547                                                                                                                                                                                                                                                                           |
| Observed Band               |                                                                                                                                                                                                                                                                                        |
| Calculated Molecular Weight |                                                                                                                                                                                                                                                                                        |
| Cell Pathway                | Nucleus .                                                                                                                                                                                                                                                                              |
| Tissue Specificity          |                                                                                                                                                                                                                                                                                        |
| Function                    | caution:The gene of ZNF547 share a non-coding exon with TRAPPC2P., function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 10 C2H2-type zinc fingers., |
| Background                  |                                                                                                                                                                                                                                                                                        |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                                                                                                                                                                                   |
| Usage suggestions           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                                                                                      |



## Products Images



Western Blot analysis of various cells using ZN547 mouse mAb