




WHAL1 mouse mAb

Catalog No	YP-mAb-12222
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WHAMMP3 WHAMML1 WHDC1L1
Immunogen	Synthesized peptide derived from human WHAL1 AA range: 95-145
Specificity	This antibody detects endogenous levels of WHAL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WHAL1;WHAMMP3 WHAMML1 WHDC1L1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using WHAL1 mouse mAb