




TBC30 mouse mAb

Catalog No	YP-mAb-12228
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TBC1D30 KIAA0984
Immunogen	Synthesized peptide derived from human TBC30 AA range: 150-200
Specificity	This antibody detects endogenous levels of TBC30 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBC30;TBC1D30 KIAA0984
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TBC30 mouse mAb