




TBCD5 mouse mAb

Catalog No	YP-mAb-12255
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TBC1D5 KIAA0210
Immunogen	Synthesized peptide derived from human TBCD5 AA range: 148-198
Specificity	This antibody detects endogenous levels of TBCD5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBCD5;TBC1D5 KIAA0210
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endosome membrane . Cytoplasmic vesicle, autophagosome . During starvation induced autophagy is relocalized from endosomal localization to LC3-positive autophagosomes. .
Tissue Specificity	
Function	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TBCD5 mouse mAb