




RLA0L mouse mAb

Catalog No	YP-mAb-12262
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPLP0P6
Immunogen	Synthesized peptide derived from human RLA0L AA range: 20-70
Specificity	This antibody detects endogenous levels of RLA0L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RLA0L;RPLP0P6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosolic large ribosomal subunit,preribosome, large subunit precursor,
Tissue Specificity	
Function	function:Ribosomal protein P0 is the functional equivalent of E.coli protein L10., similarity:Belongs to the ribosomal protein L10P family.,subunit:P0 forms a pentameric complex by interaction with dimers of P1 and P2.,subunit:P0 forms a pentameric complex by interaction with dimers of P1 and P2. Interacts with Lassa virus Z protein.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RLA0L mouse mAb