




MLIP mouse mAb

Catalog No	YP-mAb-12266
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MLIP C6orf142
Immunogen	Synthesized peptide derived from human MLIP AA range: 216-266
Specificity	This antibody detects endogenous levels of MLIP at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MLIP;MLIP C6orf142
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus envelope . Nucleus, PML body . Cell membrane, sarcolemma ; Peripheral membrane protein ; Cytoplasmic side .
Tissue Specificity	Primarily expressed in heart and skeletal muscle. Also detected in liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MLIP mouse mAb