




MRI mouse mAb

Catalog No	YP-mAb-12270
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	MRI C7orf49
Immunogen	Synthesized peptide derived from human MRI AA range: 40-90
Specificity	This antibody detects endogenous levels of MRI at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	17kD
Observed Band	
Cell Pathway	[Isoform 1]: Cytoplasm . Nucleus . Chromosome . Nuclear localization may depend upon interaction with XRCC5/Ku80 and XRCC6/Ku70 heterodimer (PubMed:24610814). Localizes to DNA damage sites (PubMed:27063109). . ; [Isoform 3]: Cytoplasm . Some nuclear localization may be due to passive diffusion. . ; [Isoform 4]: Cytoplasm . Nucleus . Nuclear localization may depend upon interaction with XRCC5/Ku80 and XRCC6/Ku70 heterodimer and increases upon etoposide treatment. .
Tissue Specificity	
Function	function:May act as a regulator of proteasome.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MRI mouse mAb