




IQCE mouse mAb

Catalog No	YP-mAb-12332
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IQCE KIAA1023
Immunogen	Synthesized peptide derived from human IQCE AA range: 50-100
Specificity	This antibody detects endogenous levels of IQCE at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IQCE;IQCE KIAA1023
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium membrane ; Peripheral membrane protein ; Cytoplasmic side . The EvC complex localizes at the base of cilia in the EvC zone of primary cilia in a EFCAB7-dependent manner. .
Tissue Specificity	
Function	similarity:Contains 2 IQ domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using IQCE mouse mAb