




AGAP9 mouse mAb

Catalog No	YP-mAb-12342
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AGAP9 CTGLF6
Immunogen	Synthesized peptide derived from human AGAP9 AA range: 4-54
Specificity	This antibody detects endogenous levels of AGAP9 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGAP9;AGAP9 CTGLF6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., function:Putative GTPase-activating protein .,miscellaneous:Encoded by one of the numerous copies of centaurin gamma-like genes clustered in the q11 region of chromosome 10.,similarity:Belongs to the centaurin gamma-like family., similarity:Contains 1 ANK repeat.,similarity:Contains 1 Arf-GAP domain., similarity:Contains 1 PH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using AGAP9 mouse mAb