




UNKL mouse mAb

Catalog No	YP-mAb-12346
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UNKL C16orf28 ZC3H5L ZC3HDC5L
Immunogen	Synthesized peptide derived from human UNKL AA range: 296-346
Specificity	This antibody detects endogenous levels of UNKL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf28;EC:2.3.2.-;KIAA0734;Putative E3 ubiquitin-protein ligase UNKL;RING finger protein unkempt-like;UNKL;UNKL C16orf28 ZC3H5L ZC3HDC5L
Observed Band	75kd
Calculated Molecular Weight	75kd
Cell Pathway	[Isoform 4]: Cytoplasm. Nucleus. Isoform 4 is primarily localized in the cytoplasm but has the ability to shuttle between the nucleus and the cytoplasm.
Tissue Specificity	
Function	similarity:Belongs to the unkempt family.,similarity:Contains 1 RING-type zinc finger.,
Background	This gene encodes a RING finger protein that may function in Rac signaling. It can bind to Brg/Brm-associated factor 60b and can promote its ubiquitination. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using UNKL mouse mAb