




CO029 mouse mAb

Catalog No	YP-mAb-12349
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C15orf29
Immunogen	Synthesized peptide derived from human CO029 AA range: 243-293
Specificity	This antibody detects endogenous levels of CO029 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CO 029;tetraspanin 8;Tetraspanin-8;TM4SF3;Tspan 8;CO029;C15orf29
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm, cytoskeleton, spindle pole . Localizes to the spindle poles only during mitosis. Sequestered to the nucleus during interphase. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CO029 mouse mAb