




TM11F mouse mAb

Catalog No	YP-mAb-12363
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMPRSS11F HATL4
Immunogen	Synthesized peptide derived from human TM11F AA range: 110-160
Specificity	This antibody detects endogenous levels of TM11F at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM11F;TMPRSS11F HATL4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	function:Probable serine protease.,similarity:Belongs to the peptidase S1 family.,similarity:Contains 1 peptidase S1 domain.,similarity:Contains 1 SEA domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TM11F mouse mAb